

Composite Astronomy Kit



MODEL - AK CAK

MAKE - AK LABS

The Composite Astronomy Kit is a comprehensive educational set designed to explore basic astronomy concepts through practical models and experiments. It helps students understand celestial bodies, the solar system, eclipses, constellations, and astronomical observations, making it ideal for STEM and physics education.

KIT SPECIFICATION

Kit Type	Composite Astronomy Kit
Subject	Astronomy & Space Science
Learning Topics	Solar System, Eclipses, Constellations & Planetary Motion
Components	Astronomy Models, Charts & Educational Accessories
Material	Durable, Educational Grade Components
Application	Astronomy, STEM & Science Education
Suitable For	Schools, Colleges, STEM Labs & Astronomy Clubs

KIT DETAILS

The **Composite Astronomy Kit** is an all-in-one learning solution featuring astronomy models and experimental tools for hands-on exploration of space science. It enables students to study the solar system, planetary motion, eclipses, constellations, and celestial phenomena through interactive activities. Designed for schools, colleges, STEM labs, and

astronomy education, it encourages scientific curiosity and experiential learning.

Included Products

	Travelling Microscope (Vernier type) in Wo High precision travelling microscope with vernier scale, ideal for school and college laboratory experiments.
	Capillary Rise Apparatus with Rising Table Precision capillary rise apparatus with adjustable rising table for measuring surface tension of liquids in physics laboratories.
	AMMETER (MO-65) Moving Coil DC with Desk stand, Range 0-1, 1.5, 3, 5, 10Amp (Any one range)
	POCKET BALANCE 0.01gm/200gm (Imported) A ultra-compact, portable electronic scale that utilizes miniature strain gauge technology to provide highly precise, on-the-go digital mass readings. Perfect for field-testing, field biology, weighing micro-reagents, and quick jewelry or component checks, this pocket-sized instrument is also widely known as a portable milligram scale, mini digital balance, or jewelry pocket scale.