

Microscope Binocular LED Light with Mechanical Stage



MODEL - BI-046/A

MAKE - AK LABS

Binocular LED Microscope is designed for clear and comfortable observation of biological specimens. Equipped with dual eyepieces, four objectives, LED illumination, and a mechanical stage, it delivers excellent optical performance for educational, clinical, and laboratory applications.

PRODUCT SPECIFICATION

Viewing Head	Binocular head with paired eyepieces
Objective Lenses	4×, 10×, 45× & 100× Oil Immersion Objectives
Eyepieces	Two Paired 10× Wide Field Eyepieces
Stage	Precision Mechanical Stage with X-Y Movement
Focusing System	Coarse & Fine Coaxial Focusing Controls
Illumination	High-Brightness Adjustable LED Light
Suitable For	Schools, Colleges, Universities, Research, Clinical & Diagnostic Laboratories

PRODUCT DETAILS

Featuring two paired 10× wide-field eyepieces, 4×, 10×, 45×, and 100× oil immersion objectives, bright LED illumination, and a precision mechanical stage, this binocular microscope ensures sharp imaging and smooth specimen positioning for schools, colleges, research, and diagnostic laboratories.