

Digital Ultrasonic Cleaner



MODEL - AK B0005

MAKE - AK LABS

Digital Ultrasonic Cleaner efficiently removes dirt, grease, and contaminants from laboratory instruments, glassware, jewelry, and precision components using high-frequency ultrasonic waves. It ensures safe, thorough, and non-abrasive cleaning for professional laboratory and industrial applications.

PRODUCT SPECIFICATION

Tank Capacity	1 Liter
Ultrasonic Frequency	28 kHz
Tank Material	SS 316 Stainless Steel
Heater Power	100 W Built-in Heater
Control Type	Analog Controls
Cleaning Method	High-frequency ultrasonic cleaning with heating
Applications	Laboratories, Medical, Dental, Research, Electronics & Precision Instrument Cleaning

PRODUCT DETAILS

Featuring a 1-liter SS 316 tank, 28 kHz ultrasonic frequency, 100 W heater, and analog controls, this cleaner provides effective cleaning performance with durable construction. It is ideal for laboratories, medical, dental, research, and quality control applications.