

Bunsen Burner with Stop Cock



MODEL - CH-016

MAKE - AK LABS

The AK Labs Bunsen Burner with Stop Cock features a 150 mm brass burner tube, 12 mm diameter, and a heavy cast iron base for maximum stability. It delivers a steady, adjustable flame, making it ideal for chemistry laboratories and educational experiments.

PRODUCT SPECIFICATION

Product Type	Bunsen Burner with Stop Cock
Brass Pipe Size	150 mm Length × 12 mm Diameter
Material	High-Quality Brass Burner with Heavy Metal Base
Gas Control	Built-in Stop Cock for Precise Flame Adjustment
Base	Heavy-duty stable base for safe operation
Applications	Heating, sterilization, and chemistry laboratory experiments
Suitable For	Schools, Colleges, Universities, Research & Industrial Laboratories

PRODUCT DETAILS

The AK Labs Heavy Base Bunsen Burner is manufactured with a precision brass burner tube (150 mm length, 12 mm diameter) and an integrated stop cock for accurate gas flow control. Its heavy-duty cast iron base provides excellent stability during heating applications. Suitable for schools, colleges, universities, and research laboratories, it

ensures safe, efficient, and consistent laboratory heating.

Key Features

- Premium brass burner tube for long service life
- Built-in stop cock for smooth gas flow adjustment
- Heavy cast iron base offers excellent stability
- Produces a clean, uniform, and adjustable flame
- Corrosion-resistant brass and chrome-plated fittings
- Compatible with standard laboratory gas tubing
- Ideal for repeated laboratory heating experiments

Applications

- Chemistry laboratory heating
- School and college practical experiments
- Sterilization of laboratory equipment
- Flame tests and combustion experiments
- Glass heating demonstrations
- Scientific research laboratories
- Educational laboratory training