

## Hot Plate 8" Round (Electric) with Energy Regulator

**MODEL** - BT-019

**MAKE** - AK LABS

Electric Hot Plate with 8-inch round heating surface is designed for uniform and controlled heating of laboratory samples. Equipped with an energy regulator, it delivers precise temperature control for chemical, biological, pharmaceutical, and educational laboratory applications.

### PRODUCT SPECIFICATION

|                     |                                                                               |
|---------------------|-------------------------------------------------------------------------------|
| Product Type        | Electric Hot Plate with Energy Regulator                                      |
| Heating Plate       | 8-inch Round Cast Iron Plate                                                  |
| Temperature Control | Adjustable Energy Regulator                                                   |
| Body Material       | Powder-Coated Heavy-Duty Metal Body                                           |
| Indicators          | Power & Heating Indicator Lamps                                               |
| Heating Performance | Fast and Uniform Heat Distribution                                            |
| Suitable For        | Schools, Colleges, Universities, Research, Chemical & Industrial Laboratories |

### PRODUCT DETAILS

Featuring an 8-inch cast iron heating plate, adjustable energy regulator, heavy-duty metal body, and indicator lamps, this electric hot plate provides fast, uniform heating with reliable performance. It is ideal for schools, colleges, research laboratories, and industrial quality control testing.