

Analytical Weight Box (Brass C.P.)

MODEL - CH-003

MAKE - AK LABS

AK Labs Analytical Weight Box (Brass C.P.) is a precision calibration weight set designed for accurate laboratory weighing and balance calibration. Supplied in a polished wooden case, it offers reliable performance for educational, scientific, and industrial measurement applications.

PRODUCT SPECIFICATION

Product Type	Analytical Weight Box (Brass C.P.)
Material	Precision Brass C.P. Weights
Capacity Options	1 mg to 100 g, 1 mg to 200 g & 1 mg to 500 g
Storage Box	Premium Wooden Box with Protective Compartments
Accessories	Forceps Included for Safe Weight Handling
Applications	Balance Calibration, Precision Weighing & Laboratory Experiments
Suitable For	Schools, Colleges, Universities, Research & Industrial Laboratories

PRODUCT DETAILS

Manufactured with precision brass C.P. weights and packed in a durable wooden storage box, this analytical weight set ensures long-term accuracy and safe storage. It is ideal for balance calibration, laboratory experiments, research work, and quality control applications.