

Dissecting Microscope - Brass Parts

MODEL - BI-059

MAKE - AK LABS

AK Labs Dissecting Microscope with superior-quality brass parts provides clear stereoscopic viewing for biological dissections and specimen observation. Supplied with both 10× and 15× eyepieces, it delivers sharp magnification for educational, laboratory, and research applications.

PRODUCT SPECIFICATION

Product Type	Dissecting Microscope
Model	BI-059/B
Construction	Superior Quality Brass Parts
Eyepieces	10× & 15× Wide Field (Interchangeable)
Stage	Glass Stage with Reversible Black & White Stage Plates
Focusing	Smooth Rack & Pinion Coarse Focusing
Suitable For	Schools, Colleges, Universities, Biology Laboratories & Research Institutes

PRODUCT DETAILS

Designed with precision brass construction, this dissecting microscope includes interchangeable 10× and 15× eyepieces, smooth focusing, and glass stage with reversible black and white stage plates. Packed in a thermocol box, it is ideal for biology practicals and specimen examination.