

Plug Key One Way with 1/2" Brass Block



MODEL - PH-260/B

MAKE - AK LABS

Plug Key One Way with 1/2-inch brass block is a reliable electrical switching device designed for making and breaking circuits during physics experiments. It ensures secure electrical connections, smooth operation, and accurate experimental results in laboratory applications.

PRODUCT SPECIFICATION	
Product Type	One Way Plug Key
Contact Material	1/2 inch Solid Brass Block
Base Material	Heavy-Duty Insulated Base
Operation	Single Pole (One Way) Manual Switching
Conductivity	Low Contact Resistance with Reliable Electrical Connection
Applications	Electrical Circuit, Physics Practical & Laboratory Experiments
Suitable For	Schools, Colleges, Universities & Physics Laboratories

PRODUCT DETAILS
Manufactured with a heavy-duty 1/2-inch solid brass block mounted on a durable insulated base, this one-way plug key offers low contact resistance, excellent conductivity, and long service life. It is ideal for schools, colleges, universities, and physics laboratories.