

Vernier Caliper



MODEL - PH-357

MAKE - AK LABS

AK Labs Vernier Caliper is a precision measuring instrument designed for accurate measurement of internal, external, depth, and step dimensions. Available in wheel type and IME type, it ensures reliable performance for laboratory, engineering, and educational applications.

PRODUCT SPECIFICATION

Product Type	Vernier Caliper
Measuring Range	0-12 cm (120 mm)
Construction Options	Chrome Plated & Stainless Steel
Available Types	Wheel Type & IME Type
Measurement Functions	Internal, External, Depth & Step Measurements
Storage	Protective Plastic Storage Case
Suitable For	Schools, Colleges, Engineering Institutes, Workshops & Laboratories

PRODUCT DETAILS

Manufactured with chrome-plated or stainless steel construction, this 12 cm Vernier Caliper offers smooth sliding movement, clear graduations, and precise measurements. Supplied in a protective storage case, it is ideal for schools, colleges, workshops, engineering institutes, and laboratories.