

Cork Borer

MODEL - PH-088

MAKE - AK LABS

AK Labs Cork Borer is a precision laboratory tool designed for cutting clean, uniform holes in corks and rubber stoppers. Manufactured from high-quality steel, it provides accurate performance for chemistry, biology, and laboratory preparation applications.

PRODUCT SPECIFICATION

Product Type	Cork Borer Set
Material	High-Quality Hardened Steel
Construction	Sharp Cutting Edges with Durable Metal Handles
Available Types	German Type, Ordinary Type & Heavy Pattern
Available Sets	Set of 3 & Set of 6
Applications	Cutting Cork & Rubber Stoppers for Laboratory Use
Suitable For	Schools, Colleges, Universities, Chemistry & Biology Laboratories

PRODUCT DETAILS

Available in German, Ordinary, and Heavy Pattern variants, this cork borer set features sharp cutting edges and durable handles for smooth operation. It is ideal for preparing laboratory stoppers used with glassware, tubing, and experimental apparatus in educational and research laboratories.