

Spirit Lamp

MODEL - CH-117

MAKE - AK LABS

The AK Labs Spirit Lamp is a reliable laboratory heating source designed for educational and scientific applications. Available in stainless steel, aluminium, and brass construction, it provides clean and controlled flame heating for chemistry experiments in schools, colleges, and laboratories.

PRODUCT SPECIFICATION

Product Name	Spirit Lamp
Brand	AK LABS
Capacity	100 ml
Material Options	Stainless Steel, Aluminium & Brass
Fuel Type	Methylated Spirit / Denatured Alcohol
Cap Type	Screw Cap with Wick Holder
Applications	Laboratory Heating & Science Practical Experiments

PRODUCT DETAILS

The AK Labs Spirit Lamp is manufactured for safe and efficient laboratory heating using methylated spirit or alcohol fuel. Featuring a leak-resistant body, screw-type cap with wick holder, and corrosion-resistant construction, it ensures dependable performance during practical experiments. Available in stainless steel, aluminium, and brass variants with a 100 ml capacity, it is ideal for physics, chemistry, biology, and educational laboratories.

Key Features

- High-quality laboratory spirit lamp
- Available in stainless steel, aluminium, and brass
- Leak-resistant construction for safe operation
- Corrosion-resistant polished finish
- Provides clean and steady flame
- Easy to refill and maintain

Applications

- Chemistry laboratory heating
- Science practical experiments
- Glassware heating
- Educational laboratory demonstrations
- Biology laboratory use
- Physics laboratory applications
- School and college laboratories