

Voltmeter (EDM-80) Square Dial, Moving Coil, DC with Desk Stand (0-3 Volt)



MODEL - PH-365/A3

MAKE - AK LABS

The AK LABS EDM-80 Square Dial Voltmeter is a precision DC moving coil instrument designed for measuring voltages up to 3 volts. Featuring a clear square dial and sturdy desk stand, it is ideal for physics laboratories and educational demonstrations.

PRODUCT SPECIFICATION	
Product Name	Ammeter (EDM-80) Square Dial
Brand	AK LABS
Type	Moving Coil DC Ammeter
Current Range	0-3 Amp DC
Dial Type	EDM-80 Square Dial
Mounting	Desk Stand
Mounting	DC Current Measurement in Laboratory Experiments

PRODUCT DETAILS
The AK LABS EDM-80 Square Dial Voltmeter is a high-quality DC moving coil measuring instrument with a 0-3 Volt range. Designed for accurate low-voltage measurements, it features an easy-to-read square dial, durable housing, and stable desk stand. Suitable for physics laboratories, schools, colleges, universities, and research institutions, it ensures reliable performance during practical experiments and electrical circuit analysis.

Key Features

- Precision DC moving coil mechanism
- Clear EDM-80 square dial for easy reading
- 0-3 Volt measurement range
- Stable desk stand for laboratory use
- Durable and impact-resistant construction
- High accuracy with smooth pointer movement
- Ideal for repeated educational and laboratory applications

Applications

- Physics Laboratory Experiments
- Electrical Circuit Testing
- Voltage Measurement
- Science Practical Demonstrations
- School & College Laboratories
- University Physics Laboratories
- Educational & Research Applications