

Jockey (Pencil Jockey)



MODEL - PH-182

MAKE - AK LABS

The AK LABS Pencil Jockey is a precision contact device designed for Meter Bridge (Wheatstone Bridge) experiments. It provides smooth electrical contact with the bridge wire, ensuring accurate resistance measurements. Its ergonomic insulated handle offers comfortable handling and dependable laboratory performance.

PRODUCT SPECIFICATION	
Product Name	Jockey (Pencil Jockey)
Brand	AK LABS
Contact Tip	Precision Brass Contact
Handle Material	Insulated High-Quality Plastic
Application	Meter Bridge / Wheatstone Bridge Experiments
Operation	Smooth Sliding Contact
Suitable For	Physics Laboratory Practical Experiments

PRODUCT DETAILS
The AK LABS Pencil Jockey is manufactured for accurate and reliable use with Meter Bridge and Wheatstone Bridge apparatus. Featuring a durable insulated handle and precision brass contact tip, it ensures smooth movement along the bridge wire without causing damage. Ideal for schools, colleges, universities, and engineering laboratories, it delivers consistent electrical contact for precise experimental results.

Key Features

- Precision brass contact tip for accurate measurements
- Ergonomic insulated handle for safe and comfortable operation
- Smooth movement along the bridge wire
- Low contact resistance for reliable readings
- Durable laboratory-grade construction
- Lightweight and easy to handle
- Ideal replacement accessory for Meter Bridge apparatus

Applications

- Meter Bridge Experiments
- Wheatstone Bridge Practicals
- Unknown Resistance Measurement
- Electrical Resistance Comparison
- Physics Laboratory Demonstrations
- School & College Science Laboratories
- Engineering Electrical Experiments